


# ZOOM

sebia

“Why you should use it”

CAPIprotect

The rinse solution for  
Capillarys and Minicap analyzers



## PRODUCT INTENDED USE

CAPIprotect is a ready to use solution, which can be used as a substitute of distilled or deionized water with Clean Protect. The product is designed to clean the capillaries of both Sebia Capillarys and Minicap instruments.



### ▶ Reference

- CAPIprotect
- PN 2061

### ▶ Availability

- All Countries
- December 2015

### ▶ Instruments

- MINICAP
- MINICAP Flex Piercing
- CAPILLARYS 2
- CAPILLARYS 2 Flex Piercing
- CAPILLARYS 2 NEONAT Fast
- CAPILLARYS 3 TERA

### ▶ Specific needs for usage

- None

## MAIN BENEFITS

- Extends the capillaries life
- Increases the MTBF (Mean Time Between Failures)
- Reduces the filter consumption
- Eliminates daily water change



**1 tank = 5L**

**Stability: 2 years**

**6.5 pH**

## CHARACTERISTICS

- Ready to use
- Stability: 2 years
- Stability on board: 1 month
- 1 pack = 2 tanks
- 1 tank = 5L